








# Surgical Caps | Technical Data Sheet

Device variants		Astron surgeon cap	Univers doctors cap	Adretta Comfort bouffant cap	Adretta Basic bouffant cap	Adretta Basic bouffant cap	
							
REF		6111500	6111501	6111502	6111503	6111504	
Design	form	astronaut	beret	bouffant	bouffant	bouffant	
	seam			hemmed elastic	chain-stitch	chain-stitch	
	colour	green	green	green	blue	green	
Cut	top side	16 x 25 cm	18.5 cm ø	oval, 49 x 53 cm	round, 53 cm ø	round, 53 cm ø	
	height/lateral	26 cm	9.5 cm				
	height/front	4 cm					
	length/ribbon	36 cm					
Material	fleece	non-woven polypropylene	non-woven polypropylene	non-woven polypropylene	non-woven polypropylene	non-woven polypropylene	
	weight/m <sup>2</sup>	28 g	28 g	20 g	14 g	14 g	
	elastic	–	–	latex-free 1145 dtex Lykra clear covered by 2 x dtex 167/30/1 polyester thread	latex-free	covered by 150 den. polyester thread	covered by 150 den. polyester thread
	sewing thread	polyester	polyester	mercifill (cotton)	polyester	polyester	
	Logistics	packaging	dispenser box with 100 pieces				
	shelf life	5 years					
	storage	protect from sunlight, dust and humidity					
Sterilization		non sterile					
Quality certificates		ISO 13485; ISO 9001					
Product classification		Medical device acc. to MDD 93/42/EEC, CE class 1					

All data refer to typical single values and may be subject to alterations.