





# Visma® | Technical Data Sheet

Surgical Face Masks Type II acc. EN 14683 (non sterile)

Device variants			Visma® earloop	Visma® tie-on
				
REF			6120600	6120601
Design			with elastic ear-loops	with tie-tapes
Dimensions	mask	length	175 +/- 5 mm	175 +/- 5 mm
		width (pleated)	95 +/- 5 mm	95 +/- 5 mm
		width (total)	180 +/- 5 mm	180 +/- 5 mm
	loop/tie tape	length	180 +/- 5 mm	370 +/- 30 mm
	nose piece	length	110 +/- 4 mm	84 +/- 4 mm
		width	2.62 +/- 0.2 mm	4 +/- 0.2 mm
		thickness	0.41 +/- 0.1 mm	0.41 +/- 0.1 mm
Logistics	packaging	dispenser box with 50 pieces		
	shelf life	5 years		
	storage	protect from sunlight, dust and humidity		
Sterilization			non sterile	
Applied standards			EN 14683:2005	
Quality certificates			ISO 13485; ISO 9001	
Product classification			Medical device acc. to MDD 93/42/EEC, CE class 1	



**Device variants**

**Visma® earloop**

**Visma® tie-on**



**Material Data**

Inner layer	material	non woven polypropylene, smooth surface	
	weight/m <sup>2</sup>	20 g	17 g
	colour	white	
Outer layer	material	non woven polypropylene	
	weight/m <sup>2</sup>	20 g	17 g
	colour	green	
Filter media	material	micro spun fibers, polypropylene	
	weight/m <sup>2</sup>	25 g	
	colour	white	
Nose piece	pliable encapsulated aluminum		
Free of	glass fibers, latex, additives		

**Functional Data**

Bacterial filtration efficiency	≥ 98%
Breathing resistance	< 29,4 Pa/cm <sup>2</sup>